




MDK Antibody

Product Code	CSB-PA003244
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21741
Immunogen	Synthesized peptide derived from the C-terminal region of Human MK.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MDK; MK1; NEGF2; Midkine; MK; Amphiregulin-associated protein; ARAP; Midgestation and kidney protein; Neurite outgrowth-promoting factor 2; Neurite outgrowth-promoting protein
Immunogen Species	Homo sapiens (Human)
Target Names	MDK

Image



Western Blot analysis of 293 cells using MK Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.