





# DUSP1 Antibody

<b>Product Code</b>	CSB-PA003245
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P28562
<b>Immunogen</b>	Synthesized peptide derived from Human MKP-1 around the non-phosphorylation site of S359.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	DUSP1; CL100; MKP1; PTPN10; VH1; Dual specificity protein phosphatase 1; Dual specificity protein phosphatase hVH1; Mitogen-activated protein kinase phosphatase 1; MAP kinase phosphatase 1; MKP-1; Protein-tyrosine phosphatase CL100
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	DUSP1

## Image



Western Blot analysis of Jurkat cells using MKP-1 Polyclonal Antibody



Western Blot analysis of various cells using  
MKP-1 Polyclonal Antibody