




BIRC7 Antibody

Product Code	CSB-PA003248
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CA5
Immunogen	Synthesized peptide derived from the Internal region of Human ML-IAP.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BIRC7; KIAP; LIVIN; MLIAP; RNF50; Baculoviral IAP repeat-containing protein 7; Kidney inhibitor of apoptosis protein; KIAP; Livin; Melanoma inhibitor of apoptosis protein; ML-IAP; RING finger protein 50
Immunogen Species	Homo sapiens (Human)
Target Names	BIRC7

Image



Western Blot analysis of COLO205 cells using ML-IAP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.