




MLK4 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA003249 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q5TCX8 |
| Immunogen | Synthesized peptide derived from the Internal region of Human MLK4. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | MLK4; KIAA1804; Mitogen-activated protein kinase kinase kinase MLK4; Mixed lineage kinase 4 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MAP3K21 |

Image



Western Blot analysis of HT29 cells using MLK4 Polyclonal Antibody