




MLK4 Antibody

Product Code	CSB-PA003249
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TCX8
Immunogen	Synthesized peptide derived from the Internal region of Human MLK4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MLK4; KIAA1804; Mitogen-activated protein kinase kinase kinase MLK4; Mixed lineage kinase 4
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K21

Image



Western Blot analysis of HT29 cells using MLK4 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.