




ZAK Antibody

Product Code	CSB-PA003251
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYL2
Immunogen	Synthesized peptide derived from the C-terminal region of Human MLTK.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MLTK; ZAK; HCCS4; Mitogen-activated protein kinase kinase kinase MLT; Human cervical cancer suppressor gene 4 protein; HCCS-4; Leucine zipper and sterile alpha motif-containing kinase; MLK-like mitogen-activated protein triple kinase; Mixe
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K20

Image



Western Blot analysis of HeLa cells using MLTK Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.