




ZAK Antibody

Product Code	CSB-PA003251
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYL2
Immunogen	Synthesized peptide derived from the C-terminal region of Human MLTK.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MLTK; ZAK; HCCS4; Mitogen-activated protein kinase kinase kinase MLT; Human cervical cancer suppressor gene 4 protein; HCCS-4; Leucine zipper-and sterile alpha motif-containing kinase; MLK-like mitogen-activated protein triple kinase; Mixe
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K20

Image



Western Blot analysis of HeLa cells using MLTK Polyclonal Antibody