




MMP19 Antibody

Product Code	CSB-PA003256
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99542
Immunogen	Synthesized peptide derived from the N-terminal region of Human MMP-19.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MMP19; MMP18; RASI; Matrix metalloproteinase-19; MMP-19; Matrix metalloproteinase RASI; Matrix metalloproteinase-18; MMP-18
Immunogen Species	Homo sapiens (Human)
Target Names	MMP19

Image



Western Blot analysis of HuvEc cells using MMP-19 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.