





MNT Antibody

Product Code	CSB-PA003261
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99583
Immunogen	Synthesized peptide derived from the Internal region of Human MNT.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MNT; BHLHD3; ROX; Max-binding protein MNT; Class D basic helix-loop-helix protein 3; bHLHd3; Myc antagonist MNT; Protein ROX
Immunogen Species	Homo sapiens (Human)
Target Names	MNT

Image



Western Blot analysis of A549 cells using MNT Polyclonal Antibody



Western Blot analysis of Jurkat cells using MNT Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.