




MOB3B Antibody

Product Code	CSB-PA003263
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86TA1
Immunogen	Synthesized peptide derived from the Internal region of Human Mob3B.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MOB3B; MOBKL2B; MOB kinase activator 3B; Mob1 homolog 2b; Mps one binder kinase activator-like 2B
Immunogen Species	Homo sapiens (Human)
Target Names	MOB3B

Image



Western Blot analysis of COLO205 cells using Mob3B Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.