




MPHOSPH9 Antibody

Product Code	CSB-PA003268
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99550
Immunogen	Synthesized peptide derived from the C-terminal region of Human MPP9.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MPHOSPH9; MPP9; M-phase phosphoprotein 9
Immunogen Species	Homo sapiens (Human)
Target Names	MPHOSPH9

Image



Western Blot analysis of A549 cells using MPP9 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.