





MRE11A Antibody

Product Code	CSB-PA003273
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49959
Immunogen	Synthesized peptide derived from Human MRE11 around the non-phosphorylation site of S264.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MRE11A; HNGS1; MRE11; Double-strand break repair protein MRE11A; Meiotic recombination 11 homolog 1; MRE11 homolog 1; Meiotic recombination 11 homolog A; MRE11 homolog A
Immunogen Species	Homo sapiens (Human)
Target Names	MRE11

Image



Western Blot analysis of COS7 cells using MRE11 Polyclonal Antibody



Western Blot analysis of K562 cells using MRE11 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.