





MRE11A Antibody

Product Code	CSB-PA003273
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49959
Immunogen	Synthesized peptide derived from Human MRE11 around the non-phosphorylation site of S264.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MRE11A; HNGS1; MRE11; Double-strand break repair protein MRE11A; Meiotic recombination 11 homolog 1; MRE11 homolog 1; Meiotic recombination 11 homolog A; MRE11 homolog A
Immunogen Species	Homo sapiens (Human)
Target Names	MRE11

Image



Western Blot analysis of COS7 cells using MRE11 Polyclonal Antibody



Western Blot analysis of K562 cells using
MRE11 Polyclonal Antibody