




MRGPRES Antibody

Product Code	CSB-PA003275
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86SM8
Immunogen	Synthesized peptide derived from the Internal region of Human MRGE.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MRGPRES; GPR167; MRGE; Mas-related G-protein coupled receptor member E; G-protein coupled receptor 167
Immunogen Species	Homo sapiens (Human)
Target Names	MRGPRES

Image



Western Blot analysis of HeLa cells using MRGE Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.