




MRGPRG Antibody

Product Code	CSB-PA003279
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86SM5
Immunogen	Synthesized peptide derived from the C-terminal region of Human MRGG.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MRGPRG; GPR169; MRGG; Mas-related G-protein coupled receptor member G; G-protein coupled receptor 169
Immunogen Species	Homo sapiens (Human)
Target Names	MRGPRG

Image



Western Blot analysis of HuvEc cells using MRGG Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.