





ABCC10 Antibody

Product Code	CSB-PA003282
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T3U5
Immunogen	Synthesized peptide derived from the Internal region of Human MRP7.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ABCC10; MRP7; SIMRP7; Multidrug resistance-associated protein 7; ATP-binding cassette sub-family C member 10
Immunogen Species	Homo sapiens (Human)
Target Names	ABCC10

Image



Western Blot analysis of COS7 cells using MRP7 Polyclonal Antibody



Western Blot analysis of Jurkat A549 cells using MRP7 Polyclonal Antibody