




MRPL21 Antibody

Product Code	CSB-PA003292
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z2W9
Immunogen	Synthesized peptide derived from the Internal region of Human MRP-L21.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MRPL21; 39S ribosomal protein L21; mitochondrial; L21mt; MRP-L21
Immunogen Species	Homo sapiens (Human)
Target Names	MRPL21

Image



Western Blot analysis of COLO205 cells using MRP-L21 Polyclonal Antibody