




MRPL46 Antibody

Product Code	CSB-PA003299
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2W6
Immunogen	Synthesized peptide derived from the Internal region of Human MRP-L46.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MRPL46; C15orf4; LIECG2; 39S ribosomal protein L46; mitochondrial; L46mt; MRP-L46; P2ECSL
Immunogen Species	Homo sapiens (Human)
Target Names	MRPL46

Image



Western Blot analysis of HeLa cells using MRP-L46 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.