





MRPS12 Antibody

Product Code	CSB-PA003305
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15235
Immunogen	Synthesized peptide derived from the Internal region of Human MRP-S12.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MRPS12; RPMS12; RPSM12; 28S ribosomal protein S12; mitochondrial; MRP-S12; S12mt; MT-RPS12
Immunogen Species	Homo sapiens (Human)
Target Names	MRPS12

Image



Western Blot analysis of SH-SY5Y cells using MRP-S12 Polyclonal Antibody



Western Blot analysis of various cells using MRP-S12 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.