





MTOR Antibody

Product Code	CSB-PA003333
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42345
Immunogen	Synthesized peptide derived from Human mTOR around the non-phosphorylation site of S2448.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MTOR; FRAP; FRAP1; FRAP2; RAFT1; RAPT1; Serine/threonine-protein kinase mTOR; FK506-binding protein 12-rapamycin complex-associated protein 1; FKBP12-rapamycin complex-associated protein; Mammalian target of rapamycin; mTOR; Mechanistic tar
Immunogen Species	Homo sapiens (Human)
Target Names	MTOR


Image



Western Blot analysis of HuvEc cells using mTOR Polyclonal Antibody



Western Blot analysis of HY926 cells using mTOR Polyclonal Antibody



Western blot analysis of customer's lysis using mTOR antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.