




MYOG Antibody

Product Code	CSB-PA003347
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15173
Immunogen	Synthesized peptide derived from the Internal region of Human Myogenin.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MYOG; BHLHC3; MYF4; Myogenin; Class C basic helix-loop-helix protein 3; bHLHc3; Myogenic factor 4; Myf-4
Immunogen Species	Homo sapiens (Human)
Target Names	MYOG

Image



Western Blot analysis of A549 cells using Myogenin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.