





# MYOM2 Antibody

<b>Product Code</b>	CSB-PA003348
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P54296
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of Human Myomesin-2.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Monkey
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	MYOM2; Myomesin-2; 165 kDa connectin-associated protein; 165 kDa titin-associated protein; M-protein; Myomesin family member 2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MYOM2

## Image



Western Blot analysis of HT29 cells using Myomesin-2 Polyclonal Antibody



Western Blot analysis of RAT-MUSCLE cells using Myomesin-2 Polyclonal Antibody

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.