




NDUFA8 Antibody

Product Code	CSB-PA003386
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51970
Immunogen	Synthesized peptide derived from the Internal region of Human NDUFA8.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NDUFA8; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 8; Complex I-19kD; CI-19kD; Complex I-PGIV; CI-PGIV; NADH-ubiquinone oxidoreductase 19 kDa subunit
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFA8

Image



Western Blot analysis of HuvEc cells using NDUFA8 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.