




NEK5 Antibody

Product Code	CSB-PA003397
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P3R8
Immunogen	Synthesized peptide derived from the Internal region of Human Nek5.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NEK5; Serine/threonine-protein kinase Nek5; Never in mitosis A-related kinase 5; NimA-related protein kinase 5
Immunogen Species	Homo sapiens (Human)
Target Names	NEK5

Image



Western Blot analysis of COLO205 cells using Nek5 Polyclonal Antibody