





IQCB1 Antibody

Product Code	CSB-PA003400
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15051
Immunogen	Synthesized peptide derived from the Internal region of Human Nephrocystin-5.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IQCB1; KIAA0036; NPHP5; OK/SW-cl.85; IQ calmodulin-binding motif-containing protein 1; Nephrocystin-5; p53 and DNA damage-regulated IQ motif protein; PIQ
Immunogen Species	Homo sapiens (Human)
Target Names	IQCB1

Image



Western Blot analysis of SH-SY5Y cells using Nephrocystin-5 Polyclonal Antibody



Western Blot analysis of 293 Jurkat cells using Nephrocystin-5 Polyclonal Antibody