





NTN1 Antibody

Product Code	CSB-PA003404
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95631
Immunogen	Synthesized peptide derived from the C-terminal region of Human Netrin-1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NTN1; NTN1L; Netrin-1; Epididymis tissue protein Li 131P
Immunogen Species	Homo sapiens (Human)
Target Names	NTN1

Image



Western Blot analysis of COS7 cells using Netrin-1 Polyclonal Antibody



Western Blot analysis of 293 K562 COS cells using Netrin-1 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.