




NFATC3 Antibody

Product Code	CSB-PA003421
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12968
Immunogen	Synthesized peptide derived from Human NFATc3 around the non-phosphorylation site of S165.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NFATC3; NFAT4; Nuclear factor of activated T-cells; cytoplasmic 3; NF-ATc3; NFATc3; NFATx; T-cell transcription factor NFAT4; NF-AT4
Immunogen Species	Homo sapiens (Human)
Target Names	NFATC3

Image



Western Blot analysis of HeLa cells using NFATC3 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.