





NFKB1 Antibody

Product Code	CSB-PA003434
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19838
Immunogen	Synthesized peptide derived from Human NFκB-p105/p50 around the non-phosphorylation site of S337.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NFKB1; Nuclear factor NF-kappa-B p105 subunit; DNA-binding factor KBF1; EBP-1; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
Immunogen Species	Homo sapiens (Human)
Target Names	NFKB1

Image



Western Blot analysis of 293 cells using NFκB-p105/p50 Polyclonal Antibody



Western Blot analysis of HELA 293T cells using NFκB-p105/p50 Polyclonal Antibody

