





NGF Antibody

Product Code	CSB-PA003446
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01138
Immunogen	Synthesized peptide derived from the Internal region of Human NGF.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NGF; NGFB; Beta-nerve growth factor; Beta-NGF
Immunogen Species	Homo sapiens (Human)
Target Names	NGF

Image



Western Blot analysis of Jurkat cells using NGF Polyclonal Antibody



Western Blot analysis of various cells using NGF Polyclonal Antibody

