





NGF Antibody

Product Code	CSB-PA003446
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01138
Immunogen	Synthesized peptide derived from the Internal region of Human NGF.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NGF; NGFB; Beta-nerve growth factor; Beta-NGF
Immunogen Species	Homo sapiens (Human)
Target Names	NGF

Image



Western Blot analysis of Jurkat cells using NGF Polyclonal Antibody



Western Blot analysis of various cells using NGF Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.