





NGFR Antibody

Product Code	CSB-PA003448
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08138
Immunogen	Synthesized peptide derived from the Internal region of Human NGFR p75.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NGFR; TNFRSF16; Tumor necrosis factor receptor superfamily member 16; Gp80-LNGFR; Low affinity neurotrophin receptor p75NTR; Low-affinity nerve growth factor receptor; NGF receptor; p75 ICD; CD antigen CD271
Immunogen Species	Homo sapiens (Human)
Target Names	NGFR

Image



Western Blot analysis of HepG2 cells using NGFR p75 Polyclonal Antibody



Western Blot analysis of K562 cells using NGFR p75 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.