





TACR3 Antibody

Product Code	CSB-PA003454
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29371
Immunogen	Synthesized peptide derived from the C-terminal region of Human NK-3R.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TACR3; NK3R; TAC3R; Neuromedin-K receptor; NKR; NK-3 receptor; NK-3R; Neurokinin B receptor; Tachykinin receptor 3
Immunogen Species	Homo sapiens (Human)
Target Names	TACR3

Image



Western Blot analysis of COLO205 cells using NK-3R Polyclonal Antibody



Western Blot analysis of A549 cells using NK-3R Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.