




NKX2-4 Antibody

Product Code	CSB-PA003456
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2Z4
Immunogen	Synthesized peptide derived from the N-terminal region of Human Nkx-2.4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NKX2-4; NKX2D; Homeobox protein Nkx-2.4; Homeobox protein NK-2 homolog D
Immunogen Species	Homo sapiens (Human)
Target Names	NKX2-4

Image



Western Blot analysis of COLO205 cells using Nkx-2.4 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.