




RTN4 Antibody

Product Code	CSB-PA003461
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQC3
Immunogen	Synthesized peptide derived from the Internal region of Human Nogo A.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RTN4; KIAA0886; NOGO; My043; SP1507; Reticulon-4; Foccen; Neurite outgrowth inhibitor; Nogo protein; Neuroendocrine-specific protein; NSP; Neuroendocrine-specific protein C homolog; RTN-x; Reticulon-5
Immunogen Species	Homo sapiens (Human)
Target Names	RTN4

Image



Western Blot analysis of HeLa cells using Nogo A Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.