





NOS2 Antibody

Product Code	CSB-PA003464
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35228
Immunogen	Synthesized peptide derived from Human NOS2 around the non-phosphorylation site of Y151.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA, IHC, IF; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NOS2

Image



Western Blot analysis of SH-SY5Y cells using NOS2 Polyclonal Antibody



Western blot analysis of 3T3 cells using NOS2 Polyclonal antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.