




NOS3 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA003466 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P29474 |
| Immunogen | Synthesized peptide derived from Human NOS3 around the non-phosphorylation site of T494. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | NOS3; Nitric oxide synthase; endothelial; Constitutive NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | NOS3 |

Image



Western Blot analysis of COLO205 cells using NOS3 Polyclonal Antibody