




NRL Antibody

Product Code	CSB-PA003487
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54845
Immunogen	Synthesized peptide derived from the N-terminal region of Human Nrl.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NRL; D14S46E; Neural retina-specific leucine zipper protein; NRL
Immunogen Species	Homo sapiens (Human)
Target Names	NRL

Image



Western Blot analysis of 293 cells using Nrl Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.