




POU3F1 Antibody

Product Code	CSB-PA003510
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03052
Immunogen	Synthesized peptide derived from the C-terminal region of Human 40822.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	POU3F1; OCT6; OTF6; POU domain; class 3, transcription factor 1; Octamer-binding protein 6; Oct-6; Octamer-binding transcription factor 6; OTF-6; POU domain transcription factor SCIP
Immunogen Species	Homo sapiens (Human)
Target Names	POU3F1

Image



Western Blot analysis of A549 cells using Oct-6 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.