




OR10D4P Antibody

Product Code	CSB-PA003514
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGN7
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 10D4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR10D4P; OR10D4; Putative olfactory receptor 10D4
Immunogen Species	Homo sapiens (Human)
Target Names	OR10D4P

Image



Western Blot analysis of COLO205 cells using Olfactory receptor 10D4 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.