




OR10G2 Antibody

Product Code	CSB-PA003515
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGC3
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 10G2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR10G2; Olfactory receptor 10G2
Immunogen Species	Homo sapiens (Human)
Target Names	OR10G2

Image



Western Blot analysis of COLO205 cells using Olfactory receptor 10G2 Polyclonal Antibody