





OR10G4 Antibody

Product Code	CSB-PA003516
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGN3
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 10G4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	OR10G4

Image



Western Blot analysis of HeLa cells using Olfactory receptor 10G4 Polyclonal Antibody



Western Blot analysis of various cells using
Olfactory receptor 10G4 Polyclonal Antibody