





OR10G4 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA003516 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8NGN3 |
| Immunogen | Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 10G4. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | OR10G4 |

Image



Western Blot analysis of HeLa cells using Olfactory receptor 10G4 Polyclonal Antibody



Western Blot analysis of various cells using Olfactory receptor 10G4 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.