




OR10G7 Antibody

Product Code	CSB-PA003518
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGN6
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 10G7.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR10G7; Olfactory receptor 10G7; Olfactory receptor OR11-283
Immunogen Species	Homo sapiens (Human)
Target Names	OR10G7

Image



Western Blot analysis of HepG2 cells using Olfactory receptor 10G7 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.