




OR11G2 Antibody

Product Code	CSB-PA003523
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGC1
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 11G2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR11G2; Olfactory receptor 11G2; Olfactory receptor OR14-34
Immunogen Species	Homo sapiens (Human)
Target Names	OR11G2

Image



Western Blot analysis of MCF7 cells using Olfactory receptor 11G2 Polyclonal Antibody