




OR1A1 Antibody

Product Code	CSB-PA003526
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P1Q5
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 1A1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR1A1; Olfactory receptor 1A1; Olfactory receptor 17-7; OR17-7; Olfactory receptor OR17-11
Immunogen Species	Homo sapiens (Human)
Target Names	OR1A1

Image



Western Blot analysis of COLO205 cells using Olfactory receptor 1A1 Polyclonal Antibody