




OR1L6 Antibody

Product Code	CSB-PA003529
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGR2
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 1L6.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR1L6; OR1L7; Olfactory receptor 1L6; Olfactory receptor 1L7; Olfactory receptor OR9-30
Immunogen Species	Homo sapiens (Human)
Target Names	OR1L6

Image



Western Blot analysis of HeLa cells using Olfactory receptor 1L6 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.