




OR2AG1/OR2AG2 Antibody

Product Code	CSB-PA003533
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H205/A6NM03
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 2AG1/2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR2AG1; OR2AG3; Olfactory receptor 2AG1; HT3; Olfactory receptor 2AG3; Olfactory receptor OR11-79; OR2AG2; OR2AG2P; Olfactory receptor 2AG2
Immunogen Species	Homo sapiens (Human)
Target Names	OR2AG1/OR2AG2

Image



Western Blot analysis of Jurkat cells using Olfactory receptor 2AG1/2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.