




OR2K2 Antibody

Product Code	CSB-PA003537
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGT1
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 2K2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR2K2; OR2AR1P; Olfactory receptor 2K2; HTPCRH06; Olfactory receptor OR9-17
Immunogen Species	Homo sapiens (Human)
Target Names	OR2K2

Image



Western Blot analysis of HT29 cells using Olfactory receptor 2K2 Polyclonal Antibody