




OR2M7 Antibody

Product Code	CSB-PA003540
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NG81
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 2M7.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR2M7; Olfactory receptor 2M7; Olfactory receptor OR1-58
Immunogen Species	Homo sapiens (Human)
Target Names	OR2M7

Image



Western Blot analysis of HepG2 cells using Olfactory receptor 2M7 Polyclonal Antibody