




OR2T2/OR2T35 Antibody

Product Code	CSB-PA003544
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IF00/Q8NGX2
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 2T2/35.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR2T2; OR2T2P; Olfactory receptor 2T2; Olfactory receptor OR1-43; OR2T35; Olfactory receptor 2T35; Olfactory receptor OR1-66
Immunogen Species	Homo sapiens (Human)
Target Names	OR2T2/OR2T35

Image



Western Blot analysis of HepG2 cells using Olfactory receptor 2T2/35 Polyclonal Antibody